

---

## Product Datasheet

**Anti-phospho-PTK9 (Tyr321)/AP [Out of stock]  
(orb893326)**

**Description**

Anti-phospho-PTK9 (Tyr321)/AP

**Species/Host**

Rabbit

**Conjugation**

AP

**Immunogen**

KLH conjugated Synthesised phosphopeptide derived from human TWF1 around the phosphorylation site of Tyr321

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

1mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

40kDa

**Uniprot ID**
[Q56JV6](#), [Q12792](#), [Q91YR1](#), [Q5RJR2](#)
**Dilution Range**

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

**Expiration Date**

12 months from date of receipt.