

Product Datasheet

**Anti-Phospho-Insulin Receptor Beta (Tyr1185)/AP
[Out of stock] (orb894741)**

Description

Anti-Phospho-Insulin Receptor

Species/Host	Rabbit
Conjugation	AP
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human INSR around the phosphorylation site of Tyr1185 [DI(p-Y)ET]
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68/152kDa
Uniprot ID	P06213, P15208, P15127
Expiration Date	12 months from date of receipt.