

Product Datasheet

Anti-phospho-EIF4EBP1(Ser111)/AP (orb895112)

Custom Reagents.

Anti-phospho-EIF4EBP1(Ser111)/AP

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	AP
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human 4EBP1 around the phosphorylation site of Ser111
Target	EIF4EBP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	13kDa
Uniprot ID	Q13541
Dilution Range	WB=
Expiration Date	12 months from date of receipt.