

Product Datasheet

Anti-Phospho-PLCG 2 (Tyr753)/PE-Cy7 [Out of stock] (orb895182)

Description

Anti-Phospho-PLCG 2 (Tyr753)/PE-Cy7

Species/Host

Rabbit

Conjugation

PE/Cy7

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human PLC gamma 2 around the phosphorylation site of Tyr753

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

139kDa

Uniprot ID
Q8CIH5, P24135, P16885
Dilution Range

IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.