

Product Datasheet

Anti-phospho-Ron(Ser1349)/PE-Cy7 (orb895540)

Description	Anti-phospho-Ron(Ser1349)/PE-Cy7
Species/Host	Rabbit
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse MST1R around the phosphorylation site of Ser1349
Target	phospho-Ron (Ser1349)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q04912
Expiration Date	12 months from date of receipt.