

Product Datasheet

phospho-Amyloid Precursor Protein (Ser730)/Cy5 antibody (orb895654)

Description

phospho-Amyloid Precursor Protein

Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Cy5
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human APP770 around the phosphorylation site of Ser730
Target	phospho-Amyloid Precursor Protein (Ser730)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	3.4-9/83kDa
Uniprot ID	P05067
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.