

Product Datasheet

**Anti-phospho-AKT1 + AKT2 + AKT3
(Thr308+Thr309+Thr305)/AP (orb895824)**

Description

Anti-phospho-AKT1 + AKT2 + AKT3

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	AP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human AKT1 around the phosphorylation site of Thr308
Target	phospho-AKT1 + AKT2 + AKT3 (Thr308+Thr309+Thr305)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31749
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.