

Product Datasheet

Anti-Histone H3 (Tri Methyl K27)/PE-Cy3 (orb898244)

Description	Anti-Histone H3 (Tri Methyl K27)/PE-Cy3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PE/Cy3
Tested Applications	FC
Immunogen	KLH conjugated synthesised methylpeptide derived from human Histone H3 around the methylation site of Tri Methyl K27
Target	Histone H3 (Tri Methyl K27)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	Histone H3.1
Uniprot ID	P68431
Expiration Date	12 months from date of receipt.