

Product Datasheet

**Anti-phospho-Tau (Thr231)/PE-Cy5.5 [Out of stock]
(orb898686)**

Description

Anti-phospho-Tau (Thr231)/PE-Cy5.5

Species/Host

Rabbit

Conjugation

PE/Cy5.5

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human Tau protein isoform 2 around the phosphorylation site of Thr231 [VR(p-T)PP]

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

52/79kDa

Uniprot ID

P10637, **P29172**, **P10636**, **P19332**

Expiration Date

12 months from date of receipt.