

Product Datasheet

**Histone H3 (acetyl K9, phospho S10) /AP antibody
(orb898806)**

Description

Histone H3 (acetyl K9, phospho S10)

Species/Host	Rabbit
Reactivity	Human
Conjugation	AP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Histone H3 around the site of Acetyl-Lys9 and phospho-Ser10
Target	Histone H3 (acetyl K9, phospho S10)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15kDa
CAS Number	Histone H3.1
Uniprot ID	P68431
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.