

Product Datasheet

Phospho-Insulin Receptor Beta (Tyr1355)/PE-Cy3 antibody (orb899307)

Description

Phospho-Insulin Receptor Beta

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy3
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Insulin Receptor alpha around the phosphorylation site of Tyr1355
Target	Phospho-Insulin Receptor Beta (Tyr1355)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68/152kDa
CAS Number	3
Uniprot ID	P06213
Expiration Date	12 months from date of receipt.