

Product Datasheet

**Anti-Phospho-p63 (Ser160+Ser162)/Cy5
(orb899377)**

Description

Anti-Phospho-p63 (Ser160+Ser162)/Cy5

Species/Host

Rabbit

Reactivity

Human

Conjugation

Cy5

Tested

FC

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human p63 around the phosphorylation site of Ser160/Ser162

Target

TP63

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

77kDa

Uniprot ID
Q9H3D4
Expiration Date

12 months from date of receipt.