

Product Datasheet

Anti-Phospho-p70 S6 Kinase Beta (Thr228)/AP (orb899867)

Description	Anti-Phospho-p70 S6 Kinase Beta (Thr228)/AP
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	AP
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human RPS6KB2 around the phosphorylation site of Thr228
Target	Phospho-p70 S6 Kinase Beta (Thr228)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UBS0
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.