

Product Datasheet

Anti-Phospho-NFkB p65(Ser468)/AP (orb899888)

Description	Anti-Phospho-NFkB p65(Ser468)/AP
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	AP
Tested Applications	ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NFkBp65 around the phosphorylation site of Ser468
Target	Phospho-NFkB p65 (Ser468)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q04206
Dilution Range	WB=1:100-200,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.