

Product Datasheet

Anti-phospho-TOP2A (Ser1106)/AP (orb900226)

Description	Anti-phospho-TOP2A (Ser1106)/AP
Species/Host	Rabbit
Reactivity	Human
Conjugation	AP
Tested Applications	ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human TOPO II Alpha around the phosphorylation site of Ser1106
Target	phospho-TOP2A (Ser1106)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11388
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.