

Product Datasheet

**Acetyl and phospho-Histone H3 (Ac-K9/p-Ser10)
antibody (FITC) [Out of stock] (orb9010)**

Description

Rabbit polyclonal to Acetyl and

Ser10) antibody conjugated to FITC

Species/Host

Rabbit

Conjugation

FITC

Tested

ICC, IF

Applications
Immunogen

KLH conjugated synthetic peptide derived from human Histone H3 around the site of Acetyl-Lys9 and phospho-Ser10

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

15kDa

Uniprot ID

Q16695, **P68431**, **Q6NXT2**, **P84243**, **P68433**, **Q71DI3**, **Q6LED0**, **P84228**

Dilution Range

ICC: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.