

Product Datasheet

Anti-phospho-ERK1 (Thr197 +Thr202)/AP [Out of stock] (orb901459)

Description

Anti-phospho-ERK1 (Thr197)

Species/Host

Rabbit

Conjugation

AP

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human ERK1 around the phosphorylation site of Thr197/Thr202

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

41kDa

Uniprot ID

P21708, P63085, P63086, Q63844, P27361, P28482

Expiration Date

12 months from date of receipt.