

Product Datasheet

phospho-VEGFR2 (Tyr1175)/PE-Cy7 antibody (orb901551)

Description	phospho-VEGFR2 (Tyr1175)/PE-Cy7 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human VEGFR2 around the phosphorylation site of Tyr1175
Target	phospho-VEGFR2 (Tyr1175)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	147kDa
Uniprot ID	P35968
Expiration Date	12 months from date of receipt.