

Product Datasheet

**Anti-Phospho-MKK3(Ser218) + MKK6(Ser207) /AP
(orb901637)**

Description

Anti-Phospho-MKK3(Ser218) +

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	AP
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MKK3 around the phosphorylation site of Ser218
Target	Phospho-MKK3(Ser218) + MKK6(Ser207)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P46734
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.