

Product Datasheet

Anti-Phospho-MARCKS (Ser152 + Ser156)/AP (orb901650)

Description

Anti-Phospho-MARCKS (Ser152 +

Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	AP
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from rat MARCKS around the phosphorylation site of Ser152/156 or Ser159/163(human)
Target	Phospho-MARCKS (Ser152 + Ser156)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P29966
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.