

Product Datasheet

Anti-Phospho-Smad3 (Ser423 + Ser425)/PE-Cy7 (orb901751)

Description	Anti-Phospho-Smad3 (Ser423 + Ser425)/PE-Cy7
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	PE/Cy7
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Smad3 around the phosphorylation site of Ser423/425
Target	Phospho-Smad3 (Ser423 + Ser425)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P84022
Expiration Date	12 months from date of receipt.