

Product Datasheet

Acetyl-Histone H4 (K12) antibody (FITC) (orb9020)

Description	Acetyl-Histone H4 (K12) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated Synthesised acetylpeptide derived from human Histone H4 around the acetylation site of Lys12
Target	Histone H4 (Acetyl K12)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62805
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.