

Product Datasheet

phospho-c-Jun (Ser243)/Cy5 antibody [Out of stock] (orb902199)

Description	phospho-c-Jun (Ser243)/Cy5 antibody
Species/Host	Rabbit
Conjugation	Cy5
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human c-Jun around the phosphorylation site of Ser243
Target	phospho-c-Jun (Ser243)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43kDa
Uniprot ID	P05412
Dilution Range	Flow-Cyt=0.2µg /Test,IF=1:100-500
Expiration Date	12 months from date of receipt.