

Product Datasheet

Anti-Phospho-Connexin 43 (Ser368)/PE-Cy3 (orb902491)

Description	Anti-Phospho-Connexin 43 (Ser368)/PE-Cy3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PE/Cy3
Tested Applications	FC, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Connexin 43 around the phosphorylation site of Ser368
Target	Phospho-Connexin 43 (Ser368)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P17302
Dilution Range	Flow-Cyt=1µg /test,ICC=1:100
Expiration Date	12 months from date of receipt.