

Product Datasheet

Anti-Phospho-PLCG 2 (Tyr1217)/PE-Cy5.5 (orb902784)

Description	Anti-Phospho-PLCG 2 (Tyr1217)/PE-Cy5.5
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PLC gamma 2 around the phosphorylation site of Tyr1217
Target	Phospho-PLCG 2 (Tyr1217)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16885
Expiration Date	12 months from date of receipt.