

## Product Datasheet

**Anti-phospho-PI3 kinase p85 alpha + gamma  
(Tyr467 + Tyr199)/PE-Cy5.5 (orb902939)**

**Description**

Anti-phospho-PI3 kinase p85

|                            |                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------|
|                            | Tyr199)/PE-Cy5.5                                                                                                            |
| <b>Species/Host</b>        | Rabbit                                                                                                                      |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                           |
| <b>Conjugation</b>         | PE/Cy5.5                                                                                                                    |
| <b>Tested Applications</b> | FC                                                                                                                          |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from human PI3 kinase p85 alpha around the phosphorylation site of Tyr467 |
| <b>Target</b>              | phospho-PI3 kinase p85 alpha + gamma (Tyr467 + Tyr199)                                                                      |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                                                       |
| <b>Form/Appearance</b>     | Liquid                                                                                                                      |
| <b>Concentration</b>       | 1mg/ml                                                                                                                      |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                                            |
| <b>Note</b>                | For research use only                                                                                                       |
| <b>Isotype</b>             | IgG                                                                                                                         |
| <b>Clonality</b>           | Polyclonal                                                                                                                  |
| <b>Uniprot ID</b>          | <b>P27986</b>                                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                             |