

Product Datasheet

Anti-phospho-ERK1/2 (Thr202 + Tyr204)/Cy5 (orb903265)

Description	Anti-phospho-ERK1/2 (Thr202 + Tyr204)/Cy5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Cy5
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse p44/42 MAPK around the phosphorylation site of Thr202/Tyr204
Target	phospho-ERK1/2 (Thr202 + Tyr204)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q63844
Dilution Range	Flow-Cyt=1µg /test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.