

Product Datasheet

Anti-phospho-c-Jun (Ser243)/PE-Cy5.5 [Out of stock] (orb904358)

Description	Anti-phospho-c-Jun (Ser243)/PE-Cy5.5
Species/Host	Rabbit
Conjugation	PE/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human c-Jun around the phosphorylation site of Ser243
Target	phospho-c-Jun (Ser243)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05412
Dilution Range	Flow-Cyt=0.2µg /Test
Expiration Date	12 months from date of receipt.