

Product Datasheet

Rabbit IgG DAB Chromogenic Kit (Blue) [Out of Stock] [Out of stock] (orb90453)

Description	Rabbit IgG DAB Chromogenic Kit (Blue)
Tested Applications	WB
Storage	-20°C for one year. DAB reagent should be protected from light.
Note	For research use only
Application notes	• Chromogenic reagent A should be stored at -20°C. If crystal appears, fully dissolve it before use. • Chromogenic working solution should be freshly prepared.
Expiration Date	Please enquire.