

Product Datasheet

Mouse IgG SABC Kit (Alkaline Phosphatase) [Out of Stock] [Out of stock] (orb90460)

Description	Mouse IgG SABC kit
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Storage	4°C for one year. Avoid freezing
Note	For research use only
Application notes	Mouse IgG refers to the animal origin of the primary antibody, not the origin of the specimen. This kit must be used on primary antibodies from mouse.
Expiration Date	12 months from date of receipt.