

Product Datasheet

Mouse IgG SABC Kit (Peroxidase) [Out of Stock] [Out of stock] (orb90472)

Description	Mouse IgG SABC kit
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Storage	4°C for one year. Avoid freezing
Note	For research use only
Application notes	Mouse IgG refers to the animal origin of the primary antibody, not the origin of the specimen. This kit must be used on primary antibodies from mouse.
Dilution Range	IHC(P): 1:50-200 IHC(F): 1:50-200 ICC: 1:50-200 Optimal dilutions should be determined by end users.
Expiration Date	12 months from date of receipt.