

Product Datasheet

Rabbit IgG SABC Kit (Peroxidase) [Out of Stock] [Out of stock] (orb90473)

Description	Rabbit IgG SABC kit
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Storage	4°C for one year. Avoid freezing
Note	For research use only
Application notes	Rabbit IgG refers to the animal origin of the primary antibody, not the origin of the specimen. This kit must be used on primary antibodies from rabbit.
Dilution Range	IHC(P): 1:50-200 IHC(F): 1:50-200 ICC: 1:50-200 Optimal dilutions should be determined by end users.
Expiration Date	12 months from date of receipt.