


Product Datasheet

Cytoplasmic and Nuclear Protein Extraction Kit (orb90513)

Description	Cytoplasmic and Nuclear Protein Extraction Kit
Storage	Upon receipt store Cytoplasmic and Nuclear Protein Extraction Kit at 4°C. It is stable at 4°C for one year. Product is shipped at ambient temperature.. Shipped in room temperature.
Note	For research use only
Expiration Date	Please enquire.



Nuclear proteins and cytoplasmic protein...



Nuclear proteins (NE) and cytoplasmic pr...