

Product Datasheet

Anti-phospho-ERK1/2 (Thr202 + Tyr204)/PE-Cy5.5 (orb905424)

Description	Anti-phospho-ERK1/2 (Thr202 + Tyr204)/PE-Cy5.5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy5.5
Tested Applications	FC, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse p44/42 MAPK around the phosphorylation site of Thr202/Tyr204
Target	phospho-ERK1/2 (Thr202 + Tyr204)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q63844
Dilution Range	Flow-Cyt=1µg /test,ICC=1:100
Expiration Date	12 months from date of receipt.