

Product Datasheet

Anti-GLUT12/PE-Cy5.5 (orb906832)

Description	Anti-GLUT12/PE-Cy5.5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy5.5
Tested Applications	ICC
Immunogen	KLH conjugated synthetic peptide derived from human GLUT12
Target	GLUT12
Preservatives	PBS, Glycerol, rAlbumin
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TD20
Dilution Range	ICC=1:100
Expiration Date	12 months from date of receipt.