

Product Datasheet

IRAK1 (phospho-Ser376) antibody (FITC) (orb9069)

Description	IRAK1 (phospho-Ser376) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IRAK1 around the phosphorylation site of Ser376
Target	Phospho-IRAK1 (Ser376)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q62406
Dilution Range	Flow-Cyt=1ug/test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.