

Product Datasheet

IRAK4 (phospho-Thr345) antibody (FITC) (orb9074)

Description

Rabbit polyclonal to IRAK4 (phospho-

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IRAK4 around the phosphorylation site of Thr345
Target	IRAK4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51kDa
Uniprot ID	Q9NWZ3
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.