

Product Datasheet

Anti-phospho-JNK1 + 2 + 3 (Thr183+Tyr185)/AP (orb907669)

Description	Anti-phospho-JNK1 + 2 + 3 (Thr183+Tyr185)/AP
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	AP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human JNK1 around the phosphorylation site of Thr183/Tyr185
Target	phospho-JNK1 + 2 + 3 (Thr183+Tyr185)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P45983
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.