

Product Datasheet

**Anti-phospho C-Myc (Thr58 / Ser62)/AP
(orb907676)**

Description

Anti-phospho C-Myc (Thr58 / Ser62)/AP

Species/Host

Rabbit

Reactivity

Rat

Conjugation

AP

Tested

IHC-P

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human c-Myc around the phosphorylation site of Thr58 and Ser62

Target

MYC

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

49kDa

Uniprot ID
P01106
Dilution Range

IHC-P=1:50-200 IHC-F=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.