

Product Datasheet

alpha Elastin/Tropoelastin/PE-Cy5.5 antibody (orb909307)

Description	alpha Elastin/Tropoelastin/PE-Cy5.5 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human Elastin
Target	Elastin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70kDa
Uniprot ID	P15502
Expiration Date	12 months from date of receipt.