

Product Datasheet

Anti-CD45RO/PE-Cy5.5 (orb909867)

Description	Anti-CD45RO/PE-Cy5.5
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE/Cy5.5
Tested Applications	FC
Immunogen	Cultured T cells from an IL-2-dependent T-cell line (CA1) prepared from human peripheral blood activated with influenza virus.
Target	CD45RO
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08575
Expiration Date	12 months from date of receipt.