

Product Datasheet

**Anti-Phospho-P38 MAPK (Thr180 + Tyr182)/AP
(orb912283)**

Description

Anti-Phospho-P38 MAPK (Thr180 +

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	AP
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human p38 MAPK around the phosphorylation site of Thr180/Tyr182
Target	Phospho-P38 MAPK (Thr180 + Tyr182)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16539
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500
Expiration Date	12 months from date of receipt.