

Product Datasheet

Anti-SALL4/AP (orb918813)

Description

Anti-SALL4/AP

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

AP

Tested

WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human SALL4 (551-650aa)

Target

SALL4

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

112kDa

Uniprot ID
[Q9UJQ4](#)
Dilution Range

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.