


Product Datasheet

YWHAE (phospho-Thr232) antibody (FITC) (orb9209)

Description	YWHAE (phospho-Thr232) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human 14-3-3E around the phosphorylation site of Thr232
Target	phospho-14-3-3 epsilon (Thr232)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62258
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
MCF7 cells
u...