

Product Datasheet

**Anti-phospho-ATF2 (Thr71 + Thr53)/AP
(orb920934)**

Description

Anti-phospho-ATF2 (Thr71 + Thr53)/AP

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

AP

Tested

ELISA, ICC, IHC-Fr, IHC-P, WB

Applications
Immunogen

KLH conjugated synthesised phosphopeptide derived from human ATF2 around the phosphorylation site of Thr71

Target

phospho-ATF2 (Thr71 + Thr53)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

55kDa

Uniprot ID

P15336

Dilution Range

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000

Expiration Date

12 months from date of receipt.