

Product Datasheet

Anti-Integrin alpha 2/CD49b/AP (orb925729)

Description	Anti-Integrin alpha 2/CD49b/AP
Species/Host	Rabbit
Reactivity	Human
Conjugation	AP
Tested Applications	ELISA, ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Integrin alpha 2/CD49b
Target	Integrin alpha 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P17301
Dilution Range	WB=1:500-2000,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.