

Product Datasheet

Anti-chemerin/AP (orb929098)

Description

Anti-chemerin/AP

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

AP

Tested

IHC-P

Applications
Immunogen

KLH conjugated synthetic peptide derived from
rabbi chemerin

Target

RARRES2

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03%
Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated
freeze/thaw cycles. The lyophilized antibody is
stable at room temperature for at least one
month and for greater than a year when kept at
-20°C. When reconstituted in sterile pH 7.4
0.01M PBS or diluent of antibody the antibody is
stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

18kDa

Uniprot ID
Q99969
Dilution Range

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200
ICC=1:50-200 not yet tested in other
applications. optimal dilutions/concentrations
should be determined by the end user.

Expiration Date

12 months from date of receipt.