

Product Datasheet

LEXSY Neo (orb93632)

Description	LEXSY Neo is an aminoglycoside antibiotic used for the selection of recombinant Leishmania strains. ...
Form/Appearance	liquid; Color: colorless
Concentration	50 mg/ml
Storage	Store at -20 °C
Note	For research use only
Application notes	Selection: For selection of recombinant Leishmania strains LEXSY Neo is added to the growth medium to a final concentration of 50 µg/ml.
Formula	C ₂₀ H ₄₀ N ₄ O ₁₀
MW	Theoretical MW: 496.6
SMILES	<chem>C[C@H]([C@@H]1[C@H]([C@@H]([C@H]([C@H](O1)O[C@@H]1[C@H](C[C@H]([C@@H]([C@H]1O)O[C@@H]1[C@@H]([C@H]([C@@](CO1)(C)O)NC)O)N)N)O)O</chem>
Expiration Date	12 months from date of receipt.