

Product Datasheet

LEXSY Tet (orb93635)

Description	Tetracycline is an antibiotic which is used at sublethal concentrations for the induction of protein...
Form/Appearance	suspension; Color: yellow
Concentration	10 mg/ml
Storage	Store at -20 °C
Note	For research use only
Formula	$C_{22}H_{24}N_2O_8$
MW	Theoretical MW: 444.4
CAS Number	64-75-5
SMILES	<chem>CC1([C@H]2C[C@H]3[C@@H](C(=O)/C(=C/N)\O)/C(=O)[C@]3(C(=O)C2=C(c2c1cccc2O)O)O)N(C)C)O</chem>
Expiration Date	12 months from date of receipt.