

Product Datasheet

Anti-phospho-14-3-3 Tau (Ser232)/Cy5 (orb938263)

Description	Anti-phospho-14-3-3 Tau (Ser232)/Cy5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Cy5
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human 14-3-3 Tau around the phosphorylation site of Ser232
Target	phospho-14-3-3 theta (Ser232)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P27348
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.